

Columbia St. Mary's

Laser Safety Test

1. Protective eyewear is always designed for a specific type of laser.

- A. True
- B. False

2. Continual low-power exposure can promote cataract formation and damage the retina.

- A. True
- B. False

3. When using a laser, it is rare that the patient's eyes must be protected.

- A. True
- B. False

4. Good fire safety includes keeping sponges and towels wet.

- A. True
- B. False

5. Lasers can be used near wet, flammable prep solutions.

- A. True
- B. False

6. An excessive oxygen level can make an airway fire more likely.

- A. True
- B. False

7. To prevent airway fires, a special laser-safe ET tube must be used.

- A. True
- B. False

8. The laser foot pedal should be placed as closely as possible to any other foot pedals.

- A. True
- B. False

9. The physician and the laser operator must communicate well.

- A. True
- B. False

Approved by MSC 02.04.15